## RosetteSep™ HLA Total Lymphocyte Enrichment Cocktail

REF	15263HLA	10 mL	For labeling 250 mL of whole blood
REF	15283HLA	40 mL	For labeling 1000 mL of whole blood

### ENGLISH

#### INTENDED USE

RosetteSep™ cell enrichment cocktails are designed for the *in vitro* enrichment of specific cell subsets from human cell sources, including whole blood.

For professional *in vitro* diagnostic use only.

#### PRODUCT SPECIFICATIONS

RosetteSep™ HLA Total Lymphocyte Enrichment Cocktail REF #15223HC.1

This cocktail contains a combination of mouse and rat monoclonal antibodies. These antibodies are bound in bispecific Tetrameric Antibody Complexes (TACs) which are directed against cell surface antigens on human hematopoietic cells (CD16, CD36, CD66b), and glycophorin A on red blood cells. The mouse monoclonal antibody subclass is IgG<sub>1</sub>. It should be kept in mind that this product is a biological reagent, and as such cannot be completely characterized or quantified. Some variability is unavoidable.

#### PRINCIPLE OF METHOD

The RosetteSep™ antibody cocktail crosslinks unwanted cells in human whole blood to multiple red blood cells (RBCs), forming immunorosettes. This increases the density of unwanted (rosetted) cells, so that they pellet along with the free RBCs when centrifuged over a buoyant density medium such as Ficoll-Paque™ PLUS (Catalog #07917/07957/07967), Ficoll-Paque™ PREMIUM (Catalog #07908/07968), or RosetteSep™ DM-L (Catalog #15705). Desired cells are never labeled with antibody and are easily collected as a highly enriched population at the interface between the plasma and the buoyant density medium.

#### QUALITY CONTROL

RosetteSep™ cell enrichment cocktails are aseptically manufactured using tightly controlled processes and pre-tested components.

Each lot of RosetteSep™ cell enrichment cocktail is sterility tested according to USP standards and Quality Control performance tested in cell separation assays using human whole blood.

#### LIMITATIONS

Product has been CE Marked in the European Union according to the 98/79/EC *In Vitro* Diagnostic Device Directive and complies with its requirements.

Additional terms and conditions may apply.

#### STORAGE

Store at 2 - 8°C. This product may be shipped at 15 - 25°C, but should be refrigerated upon receipt. Do not freeze.

#### STABILITY

Product stable at 2 - 8°C until expiry date as indicated on label.

#### PRECAUTIONS

1. For professional users only.
2. Do not use cocktail if vial contents have leaked. Unused cocktail may be disposed of according to standard laboratory procedures for non-hazardous liquids.
3. This product should be handled by trained personnel observing good laboratory practices. Once this product is added to human cells, treat the suspension as potentially biohazardous. Handling of reagents and disposal of wastes should observe all local, state, or national regulations.
4. This product is a potential irritant to eyes, respiratory system, and skin. This product may also be harmful if ingested. Avoid exposure through skin, eye contact, inhalation, and ingestion.

#### SPECIAL MATERIALS REQUIRED BUT NOT PROVIDED

##### Density Medium

Density medium refers to Ficoll-Paque™ PLUS (Catalog #07917/07957/07967), Ficoll-Paque™ PREMIUM (Catalog #07908/07968) or RosetteSep™ DM-L (Catalog #15705).

##### Recommended Medium

Dulbecco's Phosphate Buffered Saline with 2% Fetal Bovine Serum (PBS + 2% FBS, Catalog #07905).


Ficoll™, Ficoll-Paque™ PLUS and Ficoll-Paque™PREMIUM are trademarks of GE Healthcare Limited.

 **STEMCELL Technologies Inc** | 400 – 570 West 7<sup>th</sup> Avenue, Vancouver, BC | V5Z 1B3 | Canada | [www.stemcell.com](http://www.stemcell.com)

#### For Technical Assistance

Tel: +1.604.877.0713

e-mail: [info@stemcell.com](mailto:info@stemcell.com)

**MDSS GmbH**

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30175 Hannover • Germany

#29672

Version 2.1.0

March 2012

Page 1 of 15

## HANDLING AND DIRECTIONS FOR USE

Ensure that blood sample, recommended medium, density medium (see "Special Materials Required but Not Provided"), and centrifuge are all at room temperature (15 - 25°C).

1. Add RosetteSep™ HLA Cocktail at **40 µL/mL** of whole blood\* (e.g. for 2 mL of whole blood, add 80 µL of cocktail). Mix well.

*\*If using samples other than fresh whole blood, see "Notes".*

2. Incubate **20 minutes** at room temperature (15 - 25°C).
3. Dilute sample with an equal volume of PBS + 2% FBS and mix gently.
4. Layer the diluted sample on top of the density medium  
OR  
Layer the density medium underneath the diluted sample.

Be careful to minimize mixing of the density medium and sample.

See Table 1 for volume recommendations. With 50 mL conical tubes, we suggest using 15 mL of density medium, to make removal of the enriched layer after centrifugation easier.


**Table 1: Recommended Volumes and Tube Sizes**

WHOLE BLOOD (mL)	PBS + 2% FBS (mL)	DENSITY MEDIUM (mL)	TUBE SIZE (mL)
1	1	1.5	5
2	2	3	14
3	3	3	14
4	4	4	14
5	5	15	50
10	10	15	50
15	15	15	50

5. Centrifuge for **20 minutes** at 1200 x g (see "Notes") at room temperature (15 - 25°C) with the brake off.
6. Remove the enriched cells from the density medium : plasma interface.  
*Sometimes it is difficult to see the cells at the interface, especially when very rare cells are enriched. It is advisable to remove some of the density medium along with the enriched cells in order to ensure optimal recovery.*
7. Wash enriched cells once with PBS + 2% FBS. Repeat.
8. Use enriched cells as desired. We recommend that enriched samples are lysed with ammonium chloride to remove residual red blood cells prior to flow cytometric analysis (this can be done as one of the wash steps) or if residual red blood cells will interfere with subsequent assays.

## ROSETTESEP™ PROCEDURE

Numbers refer to steps in "Handling and Directions for Use".


 **STEMCELL Technologies Inc** | 400 – 570 West 7<sup>th</sup> Avenue, Vancouver, BC | V5Z 1B3 | Canada | [www.stemcell.com](http://www.stemcell.com)

### For Technical Assistance

Tel: +1.604.877.0713

e-mail: [info@stemcell.com](mailto:info@stemcell.com)

**MDSS GmbH**

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30175 Hannover • Germany

#29672

Version 2.1.0

March 2012

Page 2 of 15

## NOTES

### Samples other than Whole Blood

Although RosetteSep™ has been optimized for use with whole blood, cells can be enriched from other sources (e.g. buffy coat, leukapheresis samples). The concentration of nucleated cells in the sample should not exceed  $5 \times 10^7$  cells/mL, and red blood cells (RBCs) should be present at a ratio of at least 30 - 50 RBCs per nucleated cell.

### Conversion of g to RPM

To convert g to rpm, use the following formula:

$$RPM = \sqrt{\frac{RCF}{(1.118 \times 10^{-5}) \times (\text{Radius})}}$$

Where: RCF = relative centrifugal force (g)  
RPM = centrifuge speed in revolutions per minute  
Radius = radius of rotor in cm

### Assessing Purity

Purity of lymphocytes can be measured by flow cytometry after staining with fluorochrome conjugated antibodies against CD3 (e.g. FITC anti-CD3, Catalog #10402) and CD19 (e.g. PE anti-CD19, Catalog #10509). Lymphocytes are either CD3<sup>+</sup> or CD19<sup>+</sup> (Note that CD3<sup>+</sup>CD16<sup>+</sup> NK cells are removed during the RosetteSep™ procedure).

### Typical Results

These results are for illustrative purposes only. They were obtained using samples from normal, healthy adults. Results from individual patient samples may vary.

CATALOG #	CELL TYPE ENRICHED	PURITY
15263HLA/15283HLA	Total Lymphocytes (T (CD3 <sup>+</sup> ) & B (CD19 <sup>+</sup> ) Cells)	>80%

## TECHNICAL ASSISTANCE

For technical support please contact us by email at [info@stemcell.com](mailto:info@stemcell.com) or call +1.604.877.0713.

For more information please visit [www.stemcell.com](http://www.stemcell.com).

**EC REP MDSS GmbH**

Schiffgraben 41  
30175 Hannover, Germany


<b>REF</b> Catalog or reference number	<b>LOT</b> Batch code	 Use by: MM-YYYY
 Caution, consult accompanying documents	<b>IVD</b> In Vitro Diagnostic Medical Device	 For storage within temperature limits
<b>STERILE A</b> Aseptic Processing	 Manufacturers identification (name & address)	<b>EC REP</b> Authorized EC representative in the European Community
 CE Mark		

 **STEMCELL Technologies Inc** | 400 – 570 West 7<sup>th</sup> Avenue, Vancouver, BC | V5Z 1B3 | Canada | [www.stemcell.com](http://www.stemcell.com)

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